



Core Stabilizer Pressure Biofeedback Unit

Category: Training

Product ID: IS36

<https://www.isaansport.com/training/core-stabilizer-pressure-biofeedback-unit.php>

Pressure biofeedback unit for trunk control training and testing. Place the pad under the lower back or stomach, pump it to a starting pressure with the hand bulb, and watch the needle while you move. A steady needle means the position held; a moving needle shows the athlete and the coach straight away that it did not.



Specifications

Use	Trunk and core control exercises, training and testing
Pad	Three-section inflatable pad, 24 x 17 cm
Tube	Clear, 55 cm
Gauge	Analogue, 0 to 300 mmHg (No need to exceed 100 mmHg)
Inflation	Hand bulb with release valve
Power	None, no batteries or electronics

Manuals and Set-Up

Set-up manual (English)

<https://www.isaansport.com/manuals/EN-setup-manual-pressure-biofeedback-unit.pdf>

Set-up manual (Thai)

<https://www.isaansport.com/manuals/TH-setup-manual-pressure-biofeedback-unit.pdf>

Supplied by

Isaan Sport Technology

Tha Song Khon, Muang, Maha Sarakham 44000, Thailand

support@isaansport.com

www.isaansport.com

Email support 7 days per week

For current pricing and stock status, visit the website or email for a formal quotation - support@isaansport.com quoting the Product Name above.